



Bladen Fluorochemicals Technologies Co., Ltd

PEEK Capillary Tubing Size			
ID	Outer Diameter (mm)	Inner Diameter (mm)	Wall Thickness (mm)
1	0.4	0.3	0.05
2	0.5	0.4	0.05
3	0.79	0.25	0.27
4	0.9	0.7	0.1
5	1.2	1	0.1
6	1.25	0.9	0.175
7	1.3	0.9	0.2
8	1.35	0.95	0.2
9	1.35	1.05	0.15
10	1.4	1.2	0.1
11	1.4	1.15	0.125
12	1.5	1	0.25
13	1.5	1.35	0.075
14	1.5	1.2	0.15
15	1.55	0.93	0.31
16	1.57	1	0.285
17	1.57	0.75	0.41
18	1.57	0.5	0.535
19	1.57	0.38	0.595
20	1.57	0.25	0.66
21	1.57	0.13	0.72
22	1.57	0.1	0.735
23	1.6	0.13	0.735
24	1.6	0.25	0.675
25	1.6	0.5	0.55
26	1.6	0.75	0.425
27	1.6	1	0.3
28	1.6	1.2	0.2
29	1.75	1.4	0.175

30	1.8	1.1	0.35
31	1.8	1	0.4
32	1.8	1.5	0.15
33	1.8	1.6	0.1
34	1.9	1.7	0.1
35	1.9	1.6	0.15
36	2	1	0.5
37	2.1	1.9	0.1
38	2.1	1.9	0.1
39	2.33	1.57	0.38
40	2.6	2.4	0.1
41	2.8	2.5	0.15
42	2.8	2.4	0.2
43	3	2.67	0.165
44	3.2	1.6	0.8
45	3.2	2.1	0.55
46	3.2	1.6	0.8
47	3.2	1.2	1
48	3.3	2.9	0.2
49	4.7	2.3	1.2
50	4.95	4.15	0.4
51	4.95	4.25	0.35
52	4.95	4.55	0.2
53	5	4.2	0.4
54	5.1	1.9	1.6
55	5.9	3	1.45
56	5.9	3.9	1
57	6	4	1
58	6	3	1.5
59	6.17	5.82	0.175
60	7.2	6.9	0.15
61	7.2	6.6	0.3
62	7.2	6.3	0.45
63	7.5	4	1.75
64	7.6	7.3	0.15
65	7.6	7	0.3

66	7.6	6.7	0.45
67	8.96	8.66	0.15
68	10	4.9	2.55
69	10.2	8	1.1
70	10.22	9.82	0.2
71	10.42	9.82	0.3
72	10.56	10.26	0.15
73	10.62	9.82	0.4
74	10.82	9.82	0.5
75	10.96	10.66	0.15
76	11	8	1.5
77	11.65	10	0.825
78	12	10	1
79	12	4	4
80	12	6	3
81	12	9.5	1.25
82	12	11	0.5
83	12.2	5	3.6
84	12.45	8.88	1.785
85	13.66	12.66	0.5
86	13.66	13.06	0.3
87	17	9	4
88	18	12	3
89	18.9	6.5	6.2